

Pinpoint by Phosphorus COVID-19 Test Home Collection Kit

IN-CLINIC ORDERING GUIDE

For Prescription, In Vitro Diagnostic and Emergency Use Authorization Only

In-Clinic Test Ordering Overview

ORDERING NOTE



1 Patient visits the clinic and provider deems patient in need of testing.



2 Provider orders the test using the paper test requisition form or online practice portal. Sample is then collected by the patient in the clinical setting according to the device's included instructions.



3 Provider attaches an identifying label to the sample and places sample in biohazard bag with the patient's test requisition form.

Provider orders the test using the paper test requisition form or online practice portal. Electronic form of patient's requisition form including an electronic transfer of the test order will also be acceptable to perform the test. Sample is then collected by the patient in the clinical setting according to the device's included instructions.

The sample must be received in the laboratory within 56 hours of collection.



4 Provider ships the biohazard bags with all of the collected samples from the day to Phosphorus Diagnostics laboratory in Secaucus, NJ. Package must be shipped overnight.



5 Once ready, the ordering physician will receive a notification that the results are available. Results will be sent to the physician via encrypted email, fax or Phosphorus portal. Physician will share results with patient.

This product has not been FDA cleared or approved but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories. This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens; and, The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Follow directions for use. Keep out of reach of children. Avoid contact with skin and eyes. Wash with water if stabilizing liquid comes in contact with eyes or skin. If irritation persists, seek medical advice. Do NOT ingest stabilizing liquid. It is recommended to use eyewear when performing the collection procedure.

Please visit phosphorus.com/covid-19 to learn more about COVID-19 and what Phosphorus is doing to increase testing availability during this global pandemic.

Phosphorus



support@phosphorus.com
phosphorus.com